

# Adrenalin<sup>®</sup>

(epinephrine in sodium chloride injection)

## Ready-to-Use Premixed Bag Storage and Stability FAQ

**Q:** How should ADRENALIN<sup>®</sup> be stored?

**A:** ADRENALIN<sup>®</sup> should be stored between 20°C to 25°C (68°F to 77°F). Epinephrine is light sensitive. Protect from light and freezing. Keep in aluminum overwrap until ready to use. Discard after 24 hours of opening overwrap.

**Q:** When should the ADRENALIN<sup>®</sup> IV bag be removed from the aluminum overwrap?

**A:** The ADRENALIN<sup>®</sup> IV bag should only be removed from its aluminum overwrap at time of use. Discard after 24 hours of opening overwrap.

**Q:** Does the 24-hour stability period begin once the ADRENALIN<sup>®</sup> IV bag is removed from the aluminum overwrap or at the time it is connected to the infusion set for administration to the patient?

**A:** Do not open the aluminum overwrap until time of use. The 24-hour stability period begins immediately once the aluminum overwrap is opened.

**Q:** When should the ADRENALIN<sup>®</sup> IV bag be discarded?

**A:** The ADRENALIN<sup>®</sup> IV bag should be discarded after 24 hours of opening the aluminum overwrap. Unused portions should also be discarded after discontinuing infusion. Discard any ADRENALIN<sup>®</sup> IV bag if the aluminum overwrap has been previously opened or damaged.

**Q:** Does the ADRENALIN<sup>®</sup> IV bag need to be protected from light during administration, for example, with a light protection bag (eg, amber bag)?

**A:** Endo has not conducted studies to evaluate the safety, efficacy, and stability of ADRENALIN<sup>®</sup> with the use of a light protection bag (eg, amber bag) upon removing from the aluminum overwrap and during administration.

